TPE TUBING

When a single drop of fluid can make the difference in someone's health, you need smart fluid handling solutions that get the job done safely. Manufacturers worldwide rely on connector products that improve yield, cut costs and reduce time to market.

FlowLinX[®] TPE tubing offers a breakthrough in bioprocessing tubing solutions, combining advanced functionality with a commitment to sustainability. Engineered for seamless compatibility with AdvantaFlex[®] and C-Flex[®] tubing, FlowLinX[®] TPE ensures smooth integration into your existing operations without compromising performance. Manufactured using 40% first-generation biomass, FlowLinX[®] TPE tubing delivers a reduced carbon footprint that's 25% lower than traditional TPE tubing. Advance your process performance while contributing to a more sustainable future for bioprocessing.

REGULATORY OVERVIEW:

| Material of Construction | |
|-----------------------------------|--------------------------------|
| TPE (Thermoplastic elastomer) | |
| Biocompatibility | |
| USP <85>, Bacterial Endotoxin | |
| USP <88>, Biological Reactivity | Tests, Class VI, In Vivo |
| USP <87>, Biological Reactivity | Tests, Class VI, In Vitro |
| ISO 10993-4, Hemolysis | |
| ISO 10993-10, Irritation and Sens | sitization |
| ISO 10993-11, Systemic Toxicity | |
| Physiochemical | |
| USP <661>, Plastic Packaging S | ystems |
| EP 3.2.9, Rubber Closures for Co | ntainers |
| FDA 21CRF 177.2600, Rubber Arti | cles Intended for Repeated Use |
| Extractables | |
| Per USP <665> | |
| Particulate | |
| USP <788> | |
| Legislation | |
| REACH | |
| RoHS | |
| Conflict Mineral | |
| Shelf Life | |
| 5 Years, Non-Sterile/Non-Irradia | ited |







TPE TUBING

SUSTAINABLE INNOVATION IN EVERY CONNECTION



FlowLinX[®] TPE redefines what's possible in bioprocess tubing by blending performance excellence with environmental responsibility. Our advanced TPE tubing is manufactured with 40% first-generation biomass, offering a sustainable alternative that does not sacrifice the form, fit and function you expect from a trusted product.



By reducing its carbon footprint by 25% compared to traditional TPE, FlowLinX[®] tubing empowers companies to meet their operational goals while making measurable strides towards sustainability. For example, a small biotech using 215,000 feet (65,500 meters) of TPE tubing annually, can reduce it's emissions by 5.80 tons of CO₂ annually simply by switching to FlowLinX[®].

| Key Data | FlowLinX [®] TPE Tubing | Traditional TPE Tubing | Savings with FlowLinX® |
|--|--|----------------------------------|----------------------------|
| Material Composition | 40% First-Generation Biomass | 100% Standard Petroleum-Based | Eco-Friendly Biomass Blend |
| Carbon Footprint (per foot) | 0.077 kg CO ₂ /ft (0.254 kg CO ₂ /m) | 0.112 kg CO₂/ft (0.367 kg CO₂/m) | 25% Lower Carbon Emissions |
| Annual CO₂ Emissions (in feet) | 16.56 tons (215,000 ft) | 22.36 tons (215,000 ft) | 5.80 tons saved annually |
| Annual CO₂ Emissions (in meters) | 16.56 tons (65,231 m) | 22.36 tons (65,231 m) | 5.80 tons saved annually |
| Environmental Impact (example) ¹ | Equivalent to driving 14,758 fewer miles/9,170 fewer kilometers | Higher Emissions | Significant Reduction |

1) Environmental Protection Agency (EPA) estimates that 1 gallon of gasoline burned produces 8.89 kg (19.6 pounds) of CO₂ for a vehicle that averages 22.6 miles per gallon.

WELDING PERFORMANCE COMPARISON

FlowLinX[®] TPE tubing demonstrates reliable performance for bioprocessing applications, particularly when compared to C-Flex[®] and AdvantaFlex[®], in a series of welding studies. The study focused on welding compatibility with both gamma irradiation and autoclave sterilization methods, using a range of tubing welder parameter settings (AdvantaFlex[®] and C-Flex[®] settings). FlowLinX[®] TPE exhibited excellent welding results across these treatment types, showing consistent weld strength and integrity, making it a reliable alternative to traditional TPE options in critical bioprocessing environments. This performance, coupled with its sustainable composition and reduced carbon footprint, makes FlowLinX[®] TPE an ideal choice for the bioprocessing industry, offering both environmental benefits and reliable, high-quality performance.



TPE TUBING

WELDING PERFORMANCE COMPARISON

| Tubing Weld Combination | Welder Setting | Treatment Type ¹ | Pressure Hold Test (36 psi) ² | Average Burst Pressure (psi) ³ |
|-------------------------|----------------|-----------------------------|---|---|
| F-Flex to F-Flex | AdvantaFlex | Gamma Irradiated | Pass | 65 |
| F-Flex to F-Flex | C-Flex | Gamma Irradiated | Pass | 66 |
| F-Flex to F-Flex | AdvantaFlex | Autoclaved | Pass | 67 |
| F-Flex to F-Flex | C-Flex | Autoclaved | Pass | 72 |
| F-Flex to F-Flex | AdvantaFlex | Gamma/Autoclaved | Pass | 67 |
| F-Flex to F-Flex | C-Flex | Gamma/Autoclaved | Pass | 67 |
| F-Flex to AdvantaFlex | AdvantaFlex | Gamma Irradiated | Pass | 65 |
| F-Flex to AdvantaFlex | AdvantaFlex | Autoclaved | Pass | 70 |
| F-Flex to C-Flex | C-Flex | Gamma Irradiated | Pass | 63 |
| F-Flex to C-Flex | C-Flex | Autoclaved | Pass | 65 |
| F-Flex to AdvantaFlex | AdvantaFlex | Gamma/Autoclaved | Pass | 67 |
| F-Flex to C-Flex | C-Flex | Gamma/Autoclaved | Pass | 67 |

1) Gamma irradiation was at a dose >25 kGy and Autoclave settings were 121° C for 20 minutes

2) Test pressure was set to 36 psi (1.5x working pressure) for 2 minutes. Working pressure for 3/8" x 5/8" TPE is 24 psi.

3) Burst pressure was defined as when the tubing began to balloon, not rupture.

PART NUMBERS

| | | | | | | | COIL SIZE | |
|-----------|---------|-----------|---------|----------------|--------------|------------------------------|----------------------|-----------------------|
| ID (INCH) | ID (MM) | OD (INCH) | OD (MM) | WALL (INCH) | WALL (MM) | WORKING PRESSURE (PSI) | 50 FT (15 M) | 100 FT (30 M) |
| 1/8 | 3.2 | 1/4 | 6.4 | 1/16 | 1.6 | 30 | F-FLEX-0125-0250-C50 | F-FLEX-0125-0250-C100 |
| 3/16 | 4.8 | 5/16 | 8 | 1/16 | 1.6 | 23 | F-FLEX-0188-0313-C50 | F-FLEX-0188-0313-C100 |
| 3/16 | 4.8 | 3/8 | 9.6 | 3/32 | 2.4 | 32 | F-FLEX-0188-0375-C50 | F-FLEX-0188-0375-C100 |
| 1/4 | 6.4 | 3/8 | 9.6 | 1/16 | 1.6 | 20 | F-FLEX-0250-0375-C50 | F-FLEX-0250-0375-C100 |
| 1/4 | 6.4 | 7/16 | 11.2 | 3/32 | 2.4 | 23 | F-FLEX-0250-0438-C50 | F-FLEX-0250-0438-C100 |
| 1/4 | 6.4 | 1/2 | 12.8 | 1/8 | 3.2 | 30 | F-FLEX-0250-0500-C50 | F-FLEX-0250-0500-C100 |
| 5/16 | 8 | 7/16 | 11.2 | 1/16 | 1.6 | 21 | F-FLEX-0313-0438-C50 | F-FLEX-0313-0438-C100 |
| 5/16 | 8 | 1/2 | 12.8 | 3/32 | 2.4 | 24 | F-FLEX-0313-0500-C50 | F-FLEX-0313-0500-C100 |
| 3/8 | 9.6 | 1/2 | 12.8 | 1/16 | 1.6 | 18 | F-FLEX-0375-0500-C50 | F-FLEX-0375-0500-C100 |
| 3/8 | 9.6 | 9/16 | 14.3 | 3/32 | 2.4 | 23 | F-FLEX-0375-0563-C50 | F-FLEX-0375-0563-C100 |
| 3/8 | 9.6 | 5/8 | 15.9 | 1/8 | 3.2 | 24 | F-FLEX-0375-0625-C50 | F-FLEX-0375-0625-C100 |
| 1/2 | 12.8 | 3/4 | 19 | 1/8 | 3.2 | 22 | F-FLEX-0500-0750-C50 | F-FLEX-0500-0750-C100 |
| 5/8 | 15.9 | 7/8 | 22.3 | 1/8 | 3.2 | 18 | F-FLEX-0625-0875-C50 | F-FLEX-0625-0875-C100 |
| 3/4 | 19 | 1 | 25.4 | 1/8 | 3.2 | 17 | F-FLEX-0750-1000-C50 | F-FLEX-0750-1000-C100 |
| 3/4 | 19 | 1 1/8 | 34.9 | 3/16 | 4.8 | 18 | F-FLEX-0750-1125-C50 | F-FLEX-0750-1125-C100 |
| 1 | 25.4 | 1 3/8 | 34.9 | 3/16 | 4.8 | 18 | F-FLEX-1000-1375-C50 | F-FLEX-1000-1375-C100 |

