

DESCRIPTION:

Custom ported and optimized beta bag design for consistent product flow and high recovery. Choose from a range of hose barb and sanitary ports from 1/8" to 3". Whether you are using a DPTE® Alpha port, Biosafe® port or another port, we can provide you with a custom ready-to-use beta bag configuration that can be easily integrated into a flow path assembly for isolator transfer.



| CATALOG NUMBER | INTERNAL PART NUMBER | DESCRIPTION | PACK QTY |
|----------------|----------------------|---|----------|
| FL-GBB-0188 | 00122902 | GETINGE 190mm BAG w/ 3/16" ID OPPOSING BARB | 1 |
| FL-GBB-0250 | 00122903 | GETINGE 190mm BAG w/ 1/4" ID OPPOSING BARB | 1 |
| FL-GBB-0375 | 00122904 | GETINGE 190mm BAG w/ 3/8" ID OPPOSING BARB | 1 |
| FL-GBB-0500 | 00122905 | GETINGE 190mm BAG w/ 1/2" ID OPPOSING BARB | 1 |
| FL-GBB-MAXITC | 00122906 | GETINGE 190mm BAG w/ 1.5" TC PORT | 1 |
| FL-GBB-3TC | 00122907 | GETINGE 190mm BAG w/ 3" TC PORT | 1 |
| FL-SBB-0188 | 00122910 | SARTORIUS BIOSAFE BAG w/ 3/16" ID OPPOSING BARB | 1 |
| FL-SBB-0250 | 00122911 | SARTORIUS BIOSAFE BAG w/ 1/4" ID OPPOSING BARB | 1 |
| FL-SBB-0375 | 00122912 | SARTORIUS BIOSAFE BAG w/ 3/8" ID OPPOSING BARB | 1 |
| FL-SBB-0500 | 00122913 | SARTORIUS BIOSAFE BAG w/ 1/2" ID OPPOSING BARB | 1 |
| FL-SBB-MAXITC | 00122914 | SARTORIUS BIOSAFE BAG w/ 1.5" TC PORT | 1 |
| FL-SBB-3TC | 00122915 | SARTORIUS BIOSAFE BAG w/ 3" TC PORT | 1 |
| FL-CBB-0188 | 00122916 | CRL 190mm BAG w/ 3/16" ID OPPOSING BARB | 1 |
| FL-CBB-0250 | 00122917 | CRL 190mm BAG w/ 1/4" ID OPPOSING BARB | 1 |
| FL-CBB-0375 | 00122918 | CRL 190mm BAG w/ 3/8" ID OPPOSING BARB | 1 |
| FL-CBB-0500 | 00122920 | CRL 190mm BAG w/ 1/2" ID OPPOSING BARB | 1 |
| FL-CBB-MAXITC | 00122921 | CRL 190mm BAG w/ 1.5" TC PORT | 1 |
| FL-CBB-3TC | 00122922 | CRL 190mm BAG w/ 3" TC PORT | 1 |