# **MOLDED SANITARY ASSEMBLIES**

When a single drop of fluid can make the difference in someone's health, you need smart fluid handling solutions that get the job done safely. Manufacturers worldwide rely on connector products that improve yield, cut costs and reduce time to market.

#### TYPICAL APPLICATIONS:

- CLOSED SINGLE-USE SYSTEMS
- STERILE FILLING
- PERISTALTIC PUMPS
- SAMPLING

### DESCRIPTION:

FlowLinX molded sanitary assemblies are the optimal choice for aseptic fluid transfer, sampling, sterile filling, and other critical bioprocess applications. Made from platinum cured silicone, molded sanitary assemblies are available in mini (0.75"), maxi (1.5") and 2" sanitary versions with a variety of tubing sizes and lengths. The onepiece leak proof design features a gasket molded into the face of the sanitary which ensures proper connection alignment, reduces operator contact with the fluid path and eliminates entrapment. Choose from standard configurations or customize to your requirements.

## **SPECIFICATIONS:**

- Manufactured in an ISO 8 Cleanroom
- ISO 9001:2015 Certified Facility
- **Animal Component Free**
- Mini (0.75"), maxi (1-1/2"), and 2" sanitary sizes
- Polysulfone or stainless-steel backup cups
- Sterility assurance level 10<sup>-6</sup> per ISO 11137







# **MOLDED SANITARY ASSEMBLIES**

FA

## LIQUID SILICONE RUBBER REGULATORY OVERVIEW:

Material of Construction
Platinum Cured Silicone
Biocompatibility
USP <85>, Bacterial Endotoxin
USP <88>, Biological Reactivity Tests, Class VI, In Vivo
ISO 10993-3, Mutagenicity
ISO 10993-4, Hemolysis
ISO 10993-5, Cytotoxicity
ISO 10993-10, Irritation and Sensitization
ISO 10993-11, Systemic Toxicity

Physiochemical
USP <661>, Plastic Packaging Systems
EP 3.1.9, Plastic Packaging Systems and Materials of Construction
FDA 21 CFR 177.2600, Rubber Articles Intended for Repeated Use
Legislation
REACH, RoHS, Conflict Mineral
Extractables
Per BPOG Protocol
Shelf Life
5 Years, Non-Sterile/Non-Irradiated

6

U

### ASSEMBLY CONFIGURATOR

#### **EXAMPLE PRODUCT ID: FAS188375MMP6U**

188375

First End Type Second End Type Length **Treatment Type Tubing Type Tubing Size** Backup Cup FA: FlowLinX S: Silicone X: Maxi Sanitary X: Maxi Sanitary P: Polysulfone **6**: 6" U: non-sterile Unreinforced: Assembly Unreinforced 188375: 0.188" x 0.375" M: Mini Sanitary M: Mini Sanitary S: Stainless 12: 12" S: sterile 188438: 0.188" x 0.438" 2: 2" Sanitary 2: 2" Sanitary 24: 24" B: Silicone 0: Open End 250375: 0.250" x 0.375" **36**: 36" Braid 250438: 0.250" x 0.438" **48**: 48" Reinforced 250500: 0.250" x 0.500" 60: 60" 313438: 0.313" x 0.438" **72**: 72" 313563: 0.313" x 0.563" 84: 84" 375500: 0.375" x 0.500" 96: 96" 375563: 0.375" x 0.563" 120: 120" 375625: 0.375" x 0.625" **500750**: 0.500" x 0.750" **625875**: 0.625" x 0.875" 7501000: 0.750" x 1.000" 10001375: 1.000" x 1.375" **Braid Reinforced:** 250: 0.250" ID 375: 0.375" ID 500: 0.500" ID **625:** 0.625" ID 1000: 1.000" ID

